

<b>IDENTITY</b> <i>(As Used on Label and List)</i> <b>Interval II</b>	NOTE: Blank spaces are not permitted. If any item is not applicable, or no information is available, the space must be marked to indicate that.
--	---

**Section I**

<b>Manufacturer's Name</b> Temrex Corporation	<b>Emergency Telephone Number</b> (800) 645-1226
<b>Address</b> <i>(Number, Street, City, State, and ZIP Code)</i> 112 Albany Avenue	<b>Telephone Number for Information</b> (516) 868-6221
P.O. Box 182	<b>Date Prepared</b> January 24, 2011
Freeport, NY 11520	<b>Signature of Preparer</b> <i>(optional)</i>

**Section II – Hazardous Ingredients/Identity Information**

Hazardous Components (Specific Chemical Identity: Common Name(s)) Recommended % <i>(optional)</i>	OSHA PEL	Other Limits ACGIH TLV
<b>Components:</b> Zinc Oxide B.P.		40.0
Calcium Sulfate		30.0
Sythetic Plasticizer		15.0
Potassium Sulfate		10.0
Flavoring		Trace

**Section III – Physical/Chemical Characteristics**

<b>Boiling Point</b> N/A	<b>Specific Gravity</b> <i>(H2O = 1)</i> N/A
<b>Vapor Pressure</b> <i>(mm Hg.)</i> N/A	<b>Melting Point</b> N/A
<b>Vapor Density</b> <i>(AIR = 1)</i> N/A	<b>Evaporation Rate</b> <i>(Butyl Acetate = 1)</i> N/A
<b>Solubility in Water</b> N/A	
<b>Appearance and Odor</b>	

**Section IV – Fire and Explosion Hazard Data**

<b>Flash Point</b> <i>(Method Used)</i> N/A	<b>Flammable Limits</b> N/A	<b>LEL</b> N/A	<b>UEL</b> N/A
<b>Extinguishing Media</b> Material does not support combustion.			
<b>Special Fire Fighting Procedures</b> Non-Flammable.			
<b>Unusual Fire and Explosion Hazards</b> N/A			

**Section V – Reactivity Data**

Stability	Unstable		Conditions to Avoid
	Stable	X	
Incompatibility ( <i>Materials to Avoid</i> )			
None			
Hazardous Decomposition or Byproducts			
Negligible.			
Hazardous Polymerization	May Occur		Conditions to Avoid
	Will Not Occur	X	

**Section VI – Health Hazard Data**

Route(s) of Entry:	Inhalation?	Skin?	Ingestion?
	N/A	None observed.	N/A

Health Hazards (*Acute and Chronic*)

None

Carcinogenicity:	NTP?	IARC Monographs?	OSHA Regulated?
	N/A	N/A	N/A

Signs and Symptoms of Exposure

Slight skin/eye irritation.

Medical Conditions Generally Aggravated by Exposure

None known.

Emergency and First Aid Procedures

Contact with skin: wash promptly with soap and water.

Contact with eyes: Irrigate with clean water. If irritation persists, seek medical advice.

**Section VII – Precautions for Safe Handling and Use**

Steps to Be Taken in Case Material is Released or Spilled

Scoop up and discarded in suitable container.

Waste Disposal Method

Disposal must be carried out in accordance with all federal, state and local regulations.

Precautions to Be Taken in Handling and Storing

Store in a cool, dry area.

Other Precautions

Handling of the product exclusively by dental specialists according to instructions for use.

**Section VIII – Control Measures**

Respiratory Protection ( <i>Specify Type</i> )			
Not required.			
Ventilation	Local Exhaust	N/A	Special
	Mechanical ( <i>General</i> )	N/A	Other
Protective Gloves		Not required.	Eye Protection
			Not required if used accordingly.
Other Protective Clothing or Equipment			
Not required.			
Work/Hygienic Practices			
Avoid contact with eyes.			