

IDENTITY <i>(As Used on Label and List)</i> Temrex CR	NOTE: Blank spaces are not permitted. If any item is not applicable, or no information is available, the space must be marked to indicate that.
--	---

Section I

Manufacturer's Name Temrex Corporation	Emergency Telephone Number (800) 645-1226
Address <i>(Number, Street, City, State, and ZIP Code)</i> 112 Albany Avenue	Telephone Number for Information (516) 868-6221
P.O. Box 182	Date Prepared January 24, 2011
Freeport, NY 11520	Signature of Preparer <i>(optional)</i>

Section II – Hazardous Ingredients/Identity Information

Hazardous Components (Specific Chemical Identity: Common Name(s))	OSHA PEL	ACGIH TLV	Other Limits Recommended	% (optional)
None listed.				

Section III – Physical/Chemical Characteristics

Boiling Point N/A		Specific Gravity (H₂O = 1) Average ~ 1.5
Vapor Pressure (mm Hg.) Negligible at ambient temperatures.		Melting Point N/A
Vapor Density (AIR = 1) N/A		Evaporation Rate (Butyl Acetate = 1) N/A
Solubility in Water Virtually insoluble.		
Appearance and Odor White and brown pastes. Mild odor.		

Section IV – Fire and Explosion Hazard Data

Flash Point (Method Used) >100 degrees C	Flammable Limits N/A	LEL N/A	UEL N/A
Extinguishing Media Water, Foam, CO2			
Special Fire Fighting Procedures No specific requirement.			
Unusual Fire and Explosion Hazards None known.			

